

## Meek Clavicle Brace

Product Model : **RHA701**



**Category: Medical Equipment, Rehabilitation Supports, Fracture Aids,**

### Summary :

Looking for a reliable clavicle brace that can help you recover from a simple break or fracture? Look no further than the Meek Clavicle Brace! This product is made of high-quality materials and is designed to keep the collarbone area immobilized, allowing the bone to heal properly. Plus, our clavicle brace is exported to over 50 countries and is backed by a GWS Brand warranty. So don't wait any longer - order your Meek Clavicle Brace today!


### Specification :

**Features:** Foam with stockinette sleeve. Straps are provided with metal fasteners. Triangular pattern for two way adjustment at back.

**Applications:** Management of clavicular injuries. Acromio-clavicular and scapular separations due to sprain. Posture correction in case of stooped shoulders.

<b>Meek Clavicle Brace</b>	
<b>Sizes</b>	<b>Codes</b>
X-Small	RHA701.1
Small	RHA701.2
Medium	RHA701.3
Large	RHA701.4
X-Large	RHA701.5

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)