

## **Five Functional ICU Bed Electric, Wooden Type ABS Panel**

Product Model : **HF-102**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### **Summary :**

GWS is proud to present our Five Functional ICU Bed Electric, Wooden Type ABS Panel. This ICU bed is perfect for hospitals, emergency rooms, and other medical facilities. It features five functions: height adjustment, backrest adjustment, knee rest adjustment, Trendelenburg and reverse Trendelenburg positions, and a wooden type ABS panel. The bed is made of high-quality materials and construction that is built to last. Order your Five Functional ICU Bed Electric today!


### **Specification :**

- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 460-700H mm

Functions by Hand Remote:

- 
- (a) Hi-low (b) Backrest elevation (c) Kneerest elevation (d) Trendlenburg / Reverse Trendlenburg

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)