

## Five Functional ICU Bed Mechanical

Product Model : **HF-113**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### Summary :

If you're looking for a top-quality ICU bed, look no further than our GWS Brand Five Functional ICU Bed Mechanical (HF113). This hospital bed is perfect for use in emergency situations or in hospitals and clinics. It's made of durable materials and features five functionalities that make it ideal for intensive care patients. Plus, it's available at a price you'll love. Order yours today!

### Specification :

#### Five Functional ICU Bed Mechanical (S.S Collapsible/Detachable Side Railing)


- Two Individual Detachable/Collapsible S.S Side Railing.
- ABS Head and Foot Panel with Locking and Corner Buffer.
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.

- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 625-820H mm

Functions by Smooth Crank System: (a) Hi-low (b) Backrest elevation (c) Kneerest elevation (d) Trendlenburg / Reverse Trendlenburg

**Optional:** Four section mattress

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)