


## ICU Bed Mechanical Laminated Panels

Product Model : **HF-114**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### Summary :

GWS Brand **ICU Bed** are made with the highest quality international standards, **HF-114 ICU Bed Mechanical (Laminated Panels & S.S Side Railing)** help in recovery of a patient in the ICU and are made with much more safe and stable process, which helps motivating for patient and healthcare system.

### Specification :


#### ICU Bed Mechanical (Laminated Panels & S.S Side Railing)

- Two Individual Detachable/Collapsible S.S Side Railing.
- Stainless Steel Head & Foot Panel with Laminated Board & Corner Buffer
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors, 2 with Brakes.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.

- Pre-treated & Epoxy Powder Coated.
- Overall Size (apprx.): 2090L x 910W x 625-820H mm
- Functions by Smooth Crank System: (a) Hi-low (b) Backrest elevation (c) Kneerest elevation (d) Trendlenburg / Reverse Trendlenburg

Optional: Four Section Mattress.

*Legacy of Excellence in Health Care...*

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmed.com](mailto:info@gwsmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmed.com](http://www.gwsmed.com)