

Examination Couch with Drawer & Cabinet

Product Model : **HF151**



Category: Medical Equipment, hospital furniture, Examination Couch & Tables,


Summary :

GWS brand examination couch is perfect for any medical facility. It features a drawer and cabinet for easy storage, and it's made from durable materials that can withstand heavy use. Plus, it's comfortable and easy to operate, making it the perfect choice for any doctor or nurse. Send your Bulk Inquiry Today!

Specification :

- Made of CRCA Sheets.
- Backrest adjustable with pneumatic pump.
- Padded with foam top with high quality washable rexene.
- Provided with three drawers and three utility cabinet.
- Foldable BP apparatus tray on the side.
- Inbuilt foot step.

Legacy of Excellence in Health Care...

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 info@gwsmed.com, gwssurgicals@gmail.com  +91-124-4915548  www.gwsmed.com