

## Three Functional ICU Bed, Mechanical

Product Model : **HF115**



**Category: Medical Equipment, hospital furniture, ICU Beds,**

### Summary :


Looking for a high-quality three functional mechanical ICU bed? Look no further than GWS Surgicals LLP. Our hospital beds are available in a variety of sizes and styles to suit your needs. Plus, our medical beds are built to the highest standards for durability and reliability. Trust GWS Surgicals LLP for all your hospital bed needs!

### Specification :

- Two Individual Collapsible Side Railing.
- ABS Head & Foot Panel with Corner Buffer
- M.S Tubular Framework.
- 4 Section Top made-up of Oval perforated CRCA M.S Sheet.
- Bed Mounted on 5" Central Locking Castors, 2 with Brakes.
- S.S Telescopic IV Rod Provisions on 4 Locations.
- Urine bag holder.
- Sitting Position Adjustable.
- Pre-treated & Epoxy Powder Coated.
- Functions by Smooth Crank System: (a) Hi-low (b) Backrest elevation (c) Kneerest elevation

*Legacy of Excellence in Health Care...*

---

 Corp Office: 712, 7th Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India  
Reg Office: UG9, 1st Floor Ocus Quantum Sector 51 Gurugram 122003 (Hr.) India

 [info@gwsmmed.com](mailto:info@gwsmmed.com), [gwssurgicals@gmail.com](mailto:gwssurgicals@gmail.com)  +91-124-4915548  [www.gwsmmed.com](http://www.gwsmmed.com)